



IMMACULATE
PRECISION

Accumax Reservoirs

Accumax Reservoirs

Sophisticated procedures require manufacturing sophistication in measurement as well as storage equipments. Only this level of commitment to quality produces ground-breaking results and exceptional outcomes. Thankfully you don't have to look too far.

Accumax Reagent Reservoirs are ergonomically designed after relentless data analysis for complete ease of the user. These wide based stable reservoirs are used for multichannel pipetting. You can expect maximum retrieval and enhanced efficiency with Accumax Reservoirs.



50 ml



10 ml




25 ml

Stability & Surety



Compact Polystyrene
Body With Broad Edges For
Repetitive Filling & Firm Stand




Angled Bottom
Design For Maximum
Reagent Recovery



Pour-off Spouts To Inhibit
The Spillage While Pouring
The Reagent

Graduated Wall For
Quick Measurement Of
The Remaining Liquid



Stable & Spill Proof
Extra Large Base



Ordering Information

Accumax Reservoirs		
Cat. No.	Volume	Qty/Case
ARV-10-B	10 ml Bulk	100 x 5
ARV-10-S	10 ml sterile	40 x 5
ARV-10-S-W	10 ml Individually Wrap Sterile	50
ARV-25-B	25 ml Bulk	100
ARV-25-S	25 ml sterile	40 x 5
ARV-25-S-W	25 ml Individually Wrap Sterile	100
ARV-50-B	50 ml Bulk	100
ARV-50-S	50 ml sterile	40 x 5
ARV-50-S-W	50 ml Individually Wrap Sterile	100



Haslab GmbH – Gewerbestrasse 6 – 4105 Biel Benken
T. 031 934 54 54 F. 031 934 54 58 www.haslab.ch

Accumax

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com



ISO 9001:2015
 ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.